

CAPROLON®



CAPROLON® ist ein monofiles, langfristig resorbierbares Nahtmaterial aus Polyglycolsäure und ϵ -Caprolacton. Aufgrund des Resorptionsprofils ist das Material hervorragend für Gewebe geeignet, die eine längere Wundunterstützung benötigen. Um den ohnehin geringen Widerstand beim Gewebedurchzug noch einmal zu reduzieren ist das Material mit Resolacton beschichtet. CAPROLON® ist in violett und als ungefärbtes Nahtmaterial erhältlich.

- langfristig absorbierbar, synthetisch
- sehr hohe Reißkraft
- belastbarer und sicherer Knotensitz
- minimale Gewebereaktion




CAPROLON®



CAPROLON® is a monofilament absorbable suture based on Polyglycolacid and ϵ -Caprolacton. Due to the absorption profile the material is superior for tissue that requires long wound support. Although monofilaments have a low tissue trauma, CAPROLON® provides an additional coating. The suture is provided violet and undyed.

- *synthetic with a long period of action*
- *high tensile strength*
- *high knot security*
- *excellent biocompatibility*













In-vitro-Degradationsversuch
In-vitro trial of suture degradation

CAPROLON®

Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

Nadel / needle	Stärke/gauge USP metric		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
5/8-Kreis Rundkörper = FR <i>5/8 circle, round bodied = FR</i>						
● 2 x FR 27 	2-0	3	0,70	violett/violet	PM 2011	24
1/2-Kreis Rundkörper = HR, HRS <i>1/2 circle, round bodied = HR, HRS</i>						
● HR 12 	7-0	0.5	0,70	violett/violet	PM 20224	24
	6-0	0.7	0,70	violett/violet	PM 2021	24
	5-0	1	0,70	violett/violet	PM 2022	24
	4-0	1.5	0,70	violett/violet	PM 2023	24
● 2 x HR 12 	5-0	1	0,70	violett/violet	PM 2024	24
	4-0	1.5	0,70	violett/violet	PM 20222	24
● HRS 12 	3-0	2	0,70	violett/violet	PM 2031	24
	2-0	3	0,70	violett/violet	PM 2032	24
	0	3.5	0,70	violett/violet	PM 2033	24
● HR 17 	6-0	0.7	0,70	violett/violet	PM 2025	24
	5-0	1	0,70	violett/violet	PM 2026	24
	4-0	1.5	0,70	violett/violet	PM 2027	24
	3-0	2	0,70	violett/violet	PM 2028	24
● 2 x HR 17 	6-0	0.7	0,70	violett/violet	PM 20223	24
	5-0	1	0,70	violett/violet	PM 2029	24
	4-0	1.5	0,70	violett/violet	PM 20210	24
● HR 22 	4-0	1.5	0,70	violett/violet	PM 20211	24
	3-0	2	0,70	violett/violet	PM 20212	24
● HRS 22 	2-0	3	0,70	violett/violet	PM 2034	24
	0	3.5	0,70	violett/violet	PM 2035	24
● HR 27 	4-0	1.5	0,70	violett/violet	PM 20213	24
	3-0	2	0,70	violett/violet	PM 20214	24
	2-0	3	0,70	violett/violet	PM 20215	24
	3-0	2	6 x 0,70	violett/violet	PM 20216	12
● 2 x HR 27 					Multipack	
	4-0	1.5	0,70	violett/violet	PM 20217	24
	3-0	2	0,70	violett/violet	PM 20218	24
	2-0	3	0,70	violett/violet	PM 20219	24

CAPROLON®

Sterile Nadel-Faden-Kombination






Needle - suture - combination, sterile

Nadel / needle	Stärke/gauge		Länge/length	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HRS 27 	2-0	3	0,70	violett/violet	PM 2036	24
	0	3.5	0,70	violett/violet	PM 2037	24
	1	4	0,70	violett/violet	PM 2038	24
● HR 35 	3-0	2	0,70	violett/violet	PM 20220	24
	2-0	3	0,70	violett/violet	PM 20221	24
● HRS 35 	2-0	3	0,70	violett/violet	PM 2039	24
	0	3.5	0,70	violett/violet	PM 20310	24
	1	4	0,70	violett/violet	PM 20311	24
● HRS 40 	2-0	3	0,90	violett/violet	PM 20312	24
	0	3.5	0,90	violett/violet	PM 20313	24
	1	4	0,90	violett/violet	PM 20314	24
● HRS 40 / Schlinge / loop 	0	3.5	1,50	violett/violet	PM 20318	24
	1	4	1,50	violett/violet	PM 20319	24
● HRS 48 	0	3.5	0,90	violett/violet	PM 20315	24
	1	4	0,90	violett/violet	PM 20316	24
	2	5	0,90	violett/violet	PM 20317	24
● 2 x HRS 48 	1	4	0,90	violett/violet	PM 20321	24
● HRS 48 / Schlinge / loop 	1	4	1,80	violett/violet	PM 20320	24

CAPROLON®

Sterile Nadel-Faden-Kombination










Needle - suture - combination, sterile

Nadel / needle	Stärke/gauge USP metric		Länge/length	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
1/2-Kreis Rundkörper stumpf = HRN <i>1/2 circle, round bodied, blunt point = HRN</i>						
○ HRN 50 / Schlinge / loop	0	3.5	1,40	violett/violet	PM 2051	24
						
1/2-Kreis schneidende Rundkörper = HRT <i>1/2-circle round bodied, trocar point = HRT</i>						
▼ HRT 18	5-0	1	0,70	violett/violet	PM 2061	24
	4-0	1.5	0,70	violett/violet	PM 2062	24
	3-0	2	0,70	violett/violet	PM 2063	24
	2-0	3	0,70	violett/violet	PM 2064	24
						
▼ HRT 22	4-0	1.5	0,70	violett/violet	PM 2065	24
	3-0	2	0,70	violett/violet	PM 2066	24
						
▼ HRT 25	4-0	1.5	0,70	violett/violet	PM 2069	24
	3-0	2	0,70	violett/violet	PM 2067	24
	2-0	3	0,70	violett/violet	PM 2068	24
						
▼ HRT 50 / Schlinge / loop	1	4	1,80	violett/violet	PM 20610	24
						

CAPROLON®

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke/gauge		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
3/8-Kreis Rundkörper = DR <i>3/8 circle, round bodied = DR</i>						
● 2 x DR 10 	7-0	1.5	0,70	violett/violet	PM 2071	24
	6-0	0.7	0,70	violett/violet	PM 2072	24
● DR 12 	6-0	0.7	0,70	violett/violet	PM 2073	24
	5-0	1	0,70	violett/violet	PM 2074	24
	4-0	1.5	0,70	violett/violet	PM 2077	24
● 2 x DR 12 	6-0	0.7	0,70	violett/violet	PM 2075	24
	5-0	1	0,70	violett/violet	PM 2076	24
● DR 18 / Schlinge / loop 	4-0	1.5	0,30	violett/violet	PM 2079	24
	3-0	2	0,30	violett/violet	PM 2078	24
Tsuge-Naht / Tsuge Suture						
Gerade Rundkörper = GR <i>straight, round bodied = GR</i>						
● 2 x GR 20 	4-0	1.5	0,35	violett/violet	PM 2091	24
● 2 x GR 65 	0	3.5	0,70	violett/violet	PM 2092	24
	1	4	0,70	violett/violet	PM 2093	24
● 2 x GRW 65 	3-0	2	0,70	violett/violet	PM 2095	24
1/2-Kreis schneidende = HS 1/2-Kreis außen schneidende mit PREMIUM-Schliff = HSM <i>1/2 circle, reverse cutting = HS</i> <i>1/2 circle, cutting, with PREMIUM cut = HSM</i>						
▽ HSM 9 	5-0	1	0,45	violett/violet	PM 2121	24
▼ HS 27 	2-0	3	0,70	violett/violet	PM 2111	24
	0	3.5	0,70	violett/violet	PM 2112	24
	1	1	0,70	violett/violet	PM 2113	24

CAPROLON®

Sterile Nadel-Faden-Kombination










Needle - suture - combination, sterile

Nadel / needle	Stärke/gauge		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ HS 35 	0	3.5	0,70	violett/violet	PM 2114	24
	1	4	0,70	violett/violet	PM 2115	24
▼ HS 40 	0	3.5	0,90	violett/violet	PM 2116	24
<p>3/8-Kreis schneidende = DS, mit PREMIUM-Schliff = DSM <i>3/8 circle, reverse cutting = DS, PREMIUM cut = DSM</i></p>						
▽ DSM 11  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	violett/violet	PM 2151	24
	6-0	0.7	0,45	ungefärbt/ undyed	PM 2152	24
	5-0	1	0,45	ungefärbt/ undyed	PM 2153	24
▽ DSM 13  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	violett/violet	PM 2154	24
	4-0	1.5	0,45	violett/violet	PM 2155	24
	5-0	1	0,45	ungefärbt/ undyed	PM 2156	24
	4-0	1.5	0,45	ungefärbt/ undyed	PM 2157	24
▼ DS 18 	5-0	1	0,45	ungefärbt/ undyed	PM 2146	24
	4-0	1.5	0,45	ungefärbt/ undyed	PM 2141	24
	3-0	2	0,45	ungefärbt/ undyed	PM 2142	24
	4-0	1.5	0,70	ungefärbt/ undyed	PM 2147	24
	3-0	2	0,70	ungefärbt/ undyed	PM 2148	24
▼ DS 18 / Schlinge / loop  Tsuge-Naht / Tsuge Suture	4-0	1.5	0,70	violett/violet	PM 2149	24
▽ DSM 18  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	violett/violet	PM 21510	24
	3-0	2	0,45	violett/violet	PM 21511	24
	4-0	1.5	0,45	ungefärbt/ undyed	PM 21512	24
	3-0	2	0,45	ungefärbt/ undyed	PM 21513	24

CAPROLON®

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke/gauge		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 DSM 18 F  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	violett/violet	PM 2158	24
	4-0	1.5	0,45	ungefärbt/ undyed	PM 2159	24
 DS 24	4-0	1.5	0,45	ungefärbt/ undyed	PM 2143	24
	3-0	2	0,45	ungefärbt/ undyed	PM 2144	24
	2-0	3	0,45	ungefärbt/ undyed	PM 2145	24
 DSM 24 PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	ungefärbt/ undyed	PM 21514	24
	3-0	2	0,45	ungefärbt undyed	PM 21515	24
	2-0	3	0,45	ungefärbt/ undyed	PM 21516	24
	2-0	3	0,70	ungefärbt/ undyed	PM 21517	24
<div style="background-color: #e0e0e0; padding: 10px;"> <p>3/8-Kreis schneidende Rundkörper (Diamond-Spitze), mit PREMIUM-Schliff = DRTM <i>3/8 circle, round bodied, trocar point, PREMIUM cut = DRTM</i></p> </div>						
 DRTM 10 PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,70	violett/violet	PM 2251	24
	6-0	0.7	0,70	violett/violet	PM 2252	24
 2 x DRTM 10 PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,70	violett/violet	PM 2253	24
	0.7	6-0	0,70	violett/violet	PM 2254	24
 DRTM 12 PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,70	violett/violet	PM 2255	24
	5-0	1	0,70	violett/violet	PM 2256	24
 2 x DRTM 12 PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,70	violett/violet	PM 2257	24
	5-0	1	0,70	violett/violet	PM 2258	24
 DRT 18	4-0	1.5	0,45	violett/violet	PM 2081	24

CAPROLON®

Sterile CAPROLON® MIC-Endo Schlingen, einzeln verpackt

Sterile CAPROLON® MIC-Endo Loops, single packed

Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content
USP	metric				
3-0	2	60	violett/violet	EPM 103	12
2-0	3	60	violett/violet	EPM 104	12
0	3.5	60	violett/violet	EPM 105	12
1	4	60	violett/violet	EPM 106	12
2	5	60	violett/violet	EPM 107	12