

PGA RESORBA®



PGA RESORBA® ist ein geflochten, resorbierbares Nahtmaterial aus Polyglycolsäure. Um das Gewebetrauma zu reduzieren und den Knotenlauf zu optimieren ist das Nahtmaterial mit Resolacton beschichtet. Kombiniert mit der hohen Ausgangsreißfestigkeit ist PGA RESORBA® ein Nahtmaterial, welches den hohen Ansprüchen für den täglichen Einsatz gerecht wird. Das Produkt ist in violetter Farbe und ungefärbt erhältlich.

- mittelfristig resorbierbares Multifilament
- hohe Reißkraft
- minimales Gewebetrauma
- glatter Gewebedurchzug
- guter Knotensitz




PGA RESORBA®



*PGA Resorba® is a braided absorbable suture based on Polyglycolacid.
To reduce the tissue trauma and optimize the knot run down the suture is coated with Resolacton. In combination with the high tensile strength PGA Resorba® offers high quality for almost every kind of surgical procedure.
The product is provided in violet and undyed.*

- *midterm absorbable multifilament*
- *high tensile strength*
- *minimized tissue trauma*
- *smooth tissue passage*
- *perfect knot security*










In-vitro-Degradationsversuch
In-vitro trial of suture degradation

PGA RESORBA®

Sterile Nadel-Faden-Kombination





Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
5/8-Kreis Rundkörper = FR <i>5/8 circle, round bodied = FR</i>						
● FR 27 	3-0	2	0,70	violett/violet	PA 1011	24
	2-0	3	0,70	violett/violet	PA 1019	24
	0	3.5	0,70	violett/violet	PA 1015	24
	1	4	0,70	violett/violet	PA 10110	24
● FR 35 	0	3.5	0,70	violett/violet	PA 1013	24
	1	4	0,70	violett/violet	PA 1014	24
● 2 x FR 35 	2-0	3	0,70	violett/violet	PA 1012	24
	0	3.5	0,90	violett/violet	PA 1018	24
● FR 40 	0	3.5	0,70	violett/violet	PA 1016	24
1/2-Kreis Rundkörper = HR, HRS, HRX <i>1/2 circle, round bodied = HR, HRS, HRX</i>						
● HR 10 	5-0	1	0,70	ungefärbt/ undyed	PA 10262	24
● HR 12 	6-0	0.7	0,45	violett/violet	PA 10273	24
	5-0	1	0,45	violett/violet	PA 10274	24
	6-0	0.7	0,70	violett/violet	PA 1021	24
	5-0	1	0,70	violett/violet	PA 1022	24
	4-0	1.5	0,70	violett/violet	PA 1023	24
	6-0	0.7	0,70	ungefärbt/ undyed	PA 10242	24
	5-0	1	0,70	ungefärbt/ undyed	PA 10265	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 10266	24
● HRS 12 	3-0	2	0,70	violett/violet	PA 10328	24
	2-0	3	0,70	violett/violet	PA 10329	24
	0	3.5	0,70	violett/violet	PA 10330	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination




Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HR 17 	5-0	1	0,70	violett/violet	PA 1024	24
	4-0	1.5	0,70	violett/violet	PA 1025	24
	3-0	2	0,70	violett/violet	PA 1026	24
	2-0	3	0,70	violett/violet	PA 1027	24
	5-0	1	0,70	ungefärbt/ undyed	PA 10248	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 10243	24
	3-0	2	0,70	ungefärbt/ undyed	PA 10244	24
	2-0	3	0,70	ungefärbt/ undyed	PA 10258	24
	● HRS 17 	2-0	3	0,70	violett/violet	PA 10378
0		3.5	0,70	violett/violet	PA 10379	24
1		4	0,70	violett/violet	PA 10380	24
● HR 17 A  Abziehnadel / control release needle	5-0	1	4 x 0,45	violett/violet	PA 10256 Multipack	12
	4-0	1.5	4 x 0,45	violett/violet	PA 10287 Multipack	12
	3-0	2	4 x 0,45	violett/violet	PA 10288 Multipack	12
	4-0	1.5	8 x 0,45	violett/violet	PA 1028 Multipack	12
	3-0	2	8 x 0,45	violett/violet	PA 1029 Multipack	12
	4-0	1.5	4 x 0,70	violett/violet	PA 10284 Multipack	12
	● HR 22 	4-0	1.5	0,70	violett/violet	PA 10210
3-0		2	0,70	violett/violet	PA 10211	24
2-0		3	0,70	violett/violet	PA 10212	24
3-0		2	4 x 0,70	violett/violet	PA 10283 Multipack	12
3-0		2	0,70	ungefärbt/ undyed	PA 10254	24
2-0		3	0,70	ungefärbt/ undyed	PA 10255	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination




Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HR 22 A  Abziehnadel / control release needle	4-0	1.5	4 x 0,45	violett/violet	PA 10252 Multipack	12
	3-0	2	4 x 0,45	violett/violet	PA 10286 Multipack	12
	2-0	3	4 x 0,45	violett/violet	PA 10276 Multipack	12
	4-0	1.5	8 x 0,45	violett/violet	PA 10213 Multipack	12
	3-0	2	8 x 0,45	violett/violet	PA 10214 Multipack	12
	2-0	3	8 x 0,45	violett/violet	PA 10215 Multipack	12
	3-0	2	4 x 0,70	violett/violet	PA 10241 Multipack	12
	4-0	1.5	8 x 0,70	violett/violet	PA 10285 Multipack	12
	● HRS 22 	2-0	3	0,70	violett/violet	PA 10331
0		3.5	0,70	violett/violet	PA 10332	24
1		4	0,70	violett/violet	PA 10333	24
● HRS 22 A  Abziehnadel / control release needle	0	3.5	4 x 0,70	violett/violet	PA 10357 Multipack	12

PGA RESORBA®

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HR 27 	4-0	1.5	0,70	violett/violet	PA 10216	24
	3-0	2	0,70	violett/violet	PA 10217	24
	2-0	3	0,70	violett/violet	PA 10218	24
	0	3.5	0,70	violett/violet	PA 10219	24
	4-0	1.5	0,90	violett/violet	PA 10220	24
	3-0	2	0,90	violett/violet	PA 10221	24
	2-0	3	0,90	violett/violet	PA 10222	24
	0	3.5	0,90	violett/violet	PA 10223	24
	3-0	2	0,70	ungefärbt/ undyed	PA 10250	24
	2-0	3	0,70	ungefärbt/ undyed	PA 10251	24
	0	3.5	0,70	ungefärbt/ undyed	PA 10253	24
	2-0	3	6 x 0,70	violett/violet	PA 10249 Multipack	12
	● HR 27 A  Abziehnadel / control release needle	4-0	1.5	4 x 0,45	violett/violet	PA 10227 Multipack
3-0		2	4 x 0,45	violett/violet	PA 10228 Multipack	12
2-0		3	4 x 0,45	violett/violet	PA 10229 Multipack	12
4-0		1.5	8 x 0,45	violett/violet	PA 10224 Multipack	12
3-0		2	8 x 0,45	violett/violet	PA 10225 Multipack	12
2-0		3	8 x 0,45	violett/violet	PA 10226 Multipack	12
3-0		2	0,70	violett/violet	PA 10280	24
4-0		1.5	4 x 0,70	violett/violet	PA 10230 Multipack	12
3-0		2	4 x 0,70	violett/violet	PA 10231 Multipack	12
2-0		3	4 x 0,70	violett/violet	PA 10232 Multipack	12
0		3.5	4 x 0,70	violett/violet	PA 10261 Multipack	12
3-0		2	4 x 0,45	ungefärbt/ undyed	PA 10263 Multipack	12
2-0		3	4 x 0,45	ungefärbt/ undyed	PA 10264 Multipack	12
● HRS 27 	2-0	3	0,70	violett/violet	PA 1030	24
	0	3.5	0,70	violett/violet	PA 1031	24
	1	4	0,70	violett/violet	PA 1032	24
	2	5	0,70	violett/violet	PA 10334	24
	2-0	3	0,90	violett/violet	PA 10335	24
	0	3.5	0,90	violett/violet	PA 1033	24
	1	4	0,90	violett/violet	PA 1034	24
	2	5	0,90	violett/violet	PA 1035	24
	1	4	4 x 0,70	violett/violet	PA 10364 Multipack	12
	2	5	4 x 0,70	violett/violet	PA 10351 Multipack	12
	2-0	3	0,70	ungefärbt/ undyed	PA 10353	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination




Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HRS 27 A  Abziehnadel / control release needle	2-0	3	4 x 0,45	violett/violet	PA 10336 Multipack	12
	0	3.5	4 x 0,45	violett/violet	PA 10337 Multipack	12
	1	4	4 x 0,45	violett/violet	PA 10338 Multipack	12
	0	3.5	8 x 0,45	violett/violet	PA 10348 Multipack	12
	1	4	8 x 0,45	violett/violet	PA 10372 Multipack	12
	0	3.5	4 x 0,70	violett/violet	PA 10339 Multipack	12
	1	4	4 x 0,70	violett/violet	PA 10340 Multipack	12
	2	5	4 x 0,70	violett/violet	PA 10341 Multipack	12
	● HRX 27 	1	4	0,90	violett/violet	PA 1041
● HRX 27 A  Abziehnadel / control release needle	1	4	4 x 0,70	violett/violet	PA 1045 Multipack	12
● HR 30 	3-0	2	0,70	violett/violet	PA 10245	24
	2-0	3	0,70	violett/violet	PA 10233	24
	0	3.5	0,70	violett/violet	PA 10234	24
● HR 30 A  Abziehnadel / control release needle	0	3.5	4 x 0,70	violett/violet	PA 10270 Multipack	12
● HRX 30 	1	4	0,90	violett/violet	PA 1042	24
	2	5	0,90	violett/violet	PA 1043	24
● HR 35 	3-0	2	0,70	violett/violet	PA 10246	24
	2-0	3	0,70	violett/violet	PA 10235	24
	0	.5	0,70	violett/violet	PA 10236	24
	1	4	0,70	violett/violet	PA 10237	24
	3-0	2	0,90	violett/violet	PA 10247	24
	2-0	3	0,90	violett/violet	PA 10238	24
	0	3.5	0,90	violett/violet	PA 10239	24
	1	4	0,90	violett/violet	PA 10240	24
	3-0	2	0,70	ungefärbt/ undyed	PA 10272	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination







Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HR 35 A  Abziehnadel / control release needle	3-0	2	4 x 0,45	violett/violet	PA 10277 Multipack	12
	2-0	3	8 x 0,70	violett/violet	PA 10269 Multipack	12
	0	3.5	4 x 0,70	violett/violet	PA 10267 Multipack	12
● HRS 35 	2-0	3	0,70	violett/violet	PA 1037	24
	0	3.5	0,70	violett/violet	PA 1038	24
	1	4	0,70	violett/violet	PA 1039	24
	3-0	2	0,90	violett/violet	PA 10371	24
	2-0	3	0,90	violett/violet	PA 10310	24
	0	3.5	0,90	violett/violet	PA 10311	24
	1	4	0,90	violett/violet	PA 10312	24
	2	5	0,90	violett/violet	PA 10313	24
	2-0	3	0,90	ungefärbt/ undyed	PA 10349	24
	2	5	0,90	ungefärbt/ undyed	PA 10376	24
● HRS 35 A  Abziehnadel / control release needle	0	3.5	4 x 0,45	violett/violet	PA 10314 Multipack	12
	1	4	4 x 0,45	violett/violet	PA 10315 Multipack	12
	2	5	4 x 0,45	violett/violet	PA 10316 Multipack	12
	0	3.5	8 x 0,50	violett/violet	PA 10324 Multipack	12
	0	3.5	4 x 0,70	violett/violet	PA 10342 Multipack	12
	1	4	4 x 0,70	violett/violet	PA 10343 Multipack	12
	2	5	4 x 0,70	violett/violet	PA 10344 Multipack	12
	1	4	4 x 0,90	violett/violet	PA 10350 Multipack	12
	0	3.5	4 x 0,45	ungefärbt/ undyed	PA 10360 Multipack	12
	1	4	4 x 0,45	ungefärbt/ undyed	PA 10361 Multipack	12

PGA RESORBA®

Sterile Nadel-Faden-Kombination







Needle - suture - combination, sterile








Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HRX 35 	0	3.5	0,70	violett/violet	PA 1047	24
	1	4	0,70	violett/violet	PA 1046	24
	1	4	0,90	violett/violet	PA 1044	24
● HR40 	2-0	3	0,90	violett/violet	PA 10282	24
● HRS 40 	2-0	3	0,70	violett/violet	PA 10370	24
	0	3.5	0,70	violett/violet	PA 10317	24
	1	4	0,70	violett/violet	PA 10318	24
	0	3.5	0,90	violett/violet	PA 10319	24
	1	4	0,90	violett/violet	PA 10320	24
	2	5	0,90	violett/violet	PA 10321	24
	0	3.5	0,70	ungefärbt/ undyed	PA 10355	24
	0	3.5	0,90	ungefärbt/ undyed	PA 10362	24
	1	4	0,90	ungefärbt/ undyed	PA 10363	24
	● HRS 40 / Schlinge / loop 	1	4	1,50	violett/violet	PA 10322
● HRS 40 A  Abziehnadel / control release needle	1	1	4 x 0,45	violett/violet	PA 10323 Multipack	12
	2	5	4 x 0,45	violett/violet	PA 10345 Multipack	12
	2-0	3	8 x 0,45	violett/violet	PA 10365 Multipack	12
	1	4	8 x 0,45	violett/violet	PA 10373 Multipack	12
	2	5	4 x 0,70	violett/violet	PA 10375 Multipack	12
● HRS 48 	3-0	2	0,70	violett/violet	PA 10359	24
	2-0	3	0,70	violett/violet	PA 10368	24
	0	3.5	0,70	violett/violet	PA 10377	24
	0	3.5	0,90	violett/violet	PA 10325	24
	1	4	0,90	violett/violet	PA 10326	24
	2	5	0,90	violett/violet	PA 10327	24
	2	5	1,20	violett/violet	PA 10354	24
	1	4	0,70	ungefärbt/ undyed	PA 10356	24
	2	5	0,90	ungefärbt/ undyed	PA 10374	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile










Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HRS 48 / Schlinge / loop 	0	3.5	1,50	violett/violet	PA 10366	24
	1	4	1,50 Schlinge/loop	violett/violet	PA 10352	24
● HRS 48 A  Abziehnadel / control release needle	1	4	4 x 0,45	violett/violet	PA 10346	12
	2	5	4 x 0,45	violett/violet	PA 10347	12
	2	5	8 x 0,45	violett/violet	PA 10369	12
	1	4	4 x 0,70	violett/violet	PA 10358	12
● HR 50 	2-0	3	0,70	violett/violet	PA 10281	24
	1	4	0,70	violett/violet	PA 10278	24
	0	3.5	0,90	violett/violet	PA 10279	24
	4	6	0,90	violett/violet	PA 10271	12
● HR 70 	3-0	2	0,70	violett/violet	PA 10275	24
	2-0	3	0,70	violett/violet	PA 10257	24
	0	3.5	0,70	violett/violet	PA 10259	24
	1	4	0,70	violett/violet	PA 10260	24
	1	4	1,20	violett/violet	PA 10268	24
● HRX 75 	2	5	1,20	violett/violet	PA 1048	24
1/2-Kreis Rundkörper stumpf = HRN 1/2 circle round bodied, blunt point= HRN						
○ HRN 50 	1	4	0,90	violett/violet	PA 1051	24

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
Gerade Rundkörper = GR <i>straight, round bodied = GR</i>						
● 2 x GR 20 	4-0	1.5	0,35	violett/violet	PA 1091	24
1/2-Kreis schneidende Rundkörper (Diamond-Spitze) = HRT, HRTX <i>1/2 circle round bodied, trocar point = HRT, HRTX</i>						
▼ HRT 18 	4-0	1.5	0,70	violett/violet	PA 1060	24
	3-0	2	0,70	violett/violet	PA 1061	24
	2-0	3	0,70	violett/violet	PA 10615	24
	0	3.5	0,70	violett/violet	PA 10616	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 10621	24
	3-0	2	0,70	ungefärbt/ undyed	PA 10613	24
	2-0	3	0,70	ungefärbt/ undyed	PA 1062	24
▼ HRT 18 A <i>Abziehnadel / control release needle</i> 	2-0	3	6 x 0,45	violett/violet	PA 10620 Multipack	12
▼ HRT 22 	3-0	2	0,70	violett/violet	PA 10632	24
	2-0	3	0,70	violett/violet	PA 10633	24
	3-0	2	0,70	ungefärbt/ undyed	PA 10625	24
	2-0	3	0,70	ungefärbt/ undyed	PA 10629	24
▼ HRT 22 A <i>Abziehnadel / control release needle</i> 	2-0	3	4 x 0,45	violett/violet	PA 1068 Multipack	12
▼ HRT 25 	3-0	2	0,70	violett/violet	PA 1064	24
	2-0	3	0,70	violett/violet	PA 1065	24
	0	3.5	0,70	violett/violet	PA 1066	24
	1	4	0,90	violett/violet	PA 1067	24
	2	5	0,90	violett/violet	PA 1063	24
	3-0	2	0,70	ungefärbt/ undyed	PA 10622	24
	2-0	3	0,70	ungefärbt/ undyed	PA 10628	24
	0	3.5	4 x 0,70	violett/violet	PA 10614 Multipack	12
▼ HRT 25 A <i>Abziehnadel / control release needle</i> 	3-0	2	4 x 0,45	ungefärbt/ undyed	PA 10630 Multipack	12
	2-0	3	4 x 0,45	ungefärbt/ undyed	PA 10631 Multipack	12

PGA RESORBA®

Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 HRT 35	2-0	3	0,90	violett/violet	PA 10610	24
	0	3.5	0,90	violett/violet	PA 10611	24
	1	4	0,90	violett/violet	PA 10612	24
	2	5	0,90	violett/violet	PA 10624	24
	2-0	3	0,90	ungefärbt/ undyed	PA 10623	24
	2-0	3	1,00	ungefärbt/ undyed	PA 1069	24
 HRT 35 A Abziehnadel / control release needle	1	4	4 x 0,45	violett/violet	PA 10619 Multipack	12
 HRT 45	0	3.5	1,50 Schlinge/loop	violett/violet	PA 10626	24
	1	4	1,50 Schlinge/loop	violett/violet	PA 10627	24
 HRTX 50	1	4	0,90	violett/violet	PA 10617	24
	2	5	0,90	violett/violet	PA 10618	24
Asymptotische schneidende Rundkörper = ART <i>Asymptotic, round bodied, trocar point = ART</i>						
 ART 25	4-0	1.5	0,45	violett/violet	PA 1281	12
3/8-Kreis Rundkörper = DR, DRF <i>3/8 circle, round bodied = DR, DRF</i>						
 DR 10 F	6-0	0.7	0,45	violett/violet	PA 1071	24
 DR 12 A Abziehnadel / control release needle	3-0	2	4 x 0,45	violett/violet	PA 1072 Multipack	12
3/8-Kreis schneidende Rundkörper = DRT <i>3/8 circle, round bodied, trocar point = DRT</i>						
 DRT 12	6-0	0.7	0,45	violett/violet	PA 1081	24
	5-0	1	0,45	violett/violet	PA 1082	24
 DRT 18	4-0	1.5	0,45	ungefärbt/ undyed	PA 1083	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination






Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
3/8-Kreis Rundkörper stumpf = DRN <i>3/8 circle, round bodied, blunt point = DRN</i>						
○ DRN 64 	0	3.5	0,90	ungefärbt/ undyed	PA 1271	24
1/2-Kreis schneidende = HS, HSX mit PREMIUM-Schliff = HSM <i>1/2 circle, reverse cutting = HS, HSX PREMIUM cut = HSM</i>						
▽ HSM 13  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	ungefärbt/ undyed	PA 1120	24
▼ HS 18 	4-0	1.5	0,70	violett/violet	PA 11112	24
	3-0	2	0,70	violett/violet	PA 1113	24
	2-0	3	0,70	violett/violet	PA 1114	24
▼ HS 22 	3-0	2	0,70	violett/violet	PA 1117	24
	2-0	3	0,70	violett/violet	PA 1118	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 11118	24
	3-0	2	0,70	ungefärbt/ undyed	PA 1119	24
	2-0	3	0,70	ungefärbt/ undyed	PA 11110	24
▼ HSX 22 	2-0	3	0,70	violett/violet	PA 11117	24
▼ HS 27 	0	3.5	0,70	violett/violet	PA 11111	24
▼ HS 35 	2-0	3	0,70	violett/violet	PA 1110	24
	0	3.5	0,70	violett/violet	PA 1111	24
	1	4	0,70	violett/violet	PA 1112	24
▼ HS 40 	1	4	0,90	violett/violet	PA 11113	24
	2	5	0,90	violett/violet	PA 11114	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ HS 50 	1	4	0,70	violett/violet	PA 11116	24
	4	6	0,90	violett/violet	PA 11115	12
▼ HS 120 	2	5	0,90	violett/violet	PA 1115	24
<div style="background-color: #cccccc; padding: 5px;"> 3/8-Kreis schneidende = DS mit PREMIUM-Schliff = DSM <i>3/8 circle, reverse cutting = DS PREMIUM cut = DSM</i> </div>						
▽ DSM 11  PREMIUM-Schliff PREMIUM cut	6-0	0.7	0,45	ungefärbt/ undyed	PA 1158	24
▼ DS 13 	4-0	1.5	0,70	ungefärbt/ undyed	PA 11422	24
▽ DSM 13  PREMIUM-Schliff PREMIUM cut	5-0	1	0,45	violett/violet	PA 1159	24
	6-0	0.7	0,70	violett/violet	PA 11512	24
	5-0	1	0,70	violett/violet	PA 11511	24
	6-0	0.7	0,45	ungefärbt/ undyed	PA 1150	24
	5-0	1	0,45	ungefärbt/ undyed	PA 1151	24
	4-0	1.5	0,45	ungefärbt/ undyed	PA 1152	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

















Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▽ 2 x DSM 13  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	violett/violet	PA 11515	24
	4-0	1.5	0,45	violett/violet	PA 11516	24
▽ 2 x DSM 13 S  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	violett/violet	PA 11514	24
▽ DSM 16  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	ungefärbt/ undyed	PA 11518	24
▼ DS 18 F 	4-0	1.5	0,45	ungefärbt/ undyed	PA 11415	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 11414	24
▼ DS 18 	5-0	1	0,45	violett/violet	PA 11417	24
	4-0	1.5	0,45	violett/violet	PA 11418	24
	4-0	1.5	0,70	violett/violet	PA 11423	24
	3-0	2	0,70	violett/violet	PA 11416	24
	5-0	1	0,45	ungefärbt/ undyed	PA 1140	24
	4-0	1.5	0,45	ungefärbt/ undyed	PA 1141	24
	3-0	2	0,45	ungefärbt/ undyed	PA 1142	24
	3-0	2	0,70	ungefärbt/ undyed	PA 11413	24
▽ DSM 18  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,70	violett/violet	PA 11510	24
	4-0	1.5	0,45	ungefärbt/ undyed	PA 1154	24
	3-0	2	0,45	ungefärbt/ undyed	PA 1155	24
▽ DSM 18 F  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	ungefärbt/ undyed	PA 1153	24

PGA RESORBA®

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ DS 24 	2-0	3	0,70	violett/violet	PA 11412	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 1143	24
	3-0	2	0,70	ungefärbt/ undyed	PA 1145	24
	2-0	3	0,70	ungefärbt/ undyed	PA 1146	24
▽ DSM 24  PREMIUM-Schliff PREMIUM cut	3-0	2	0,45	violett/violet	PA 11513	24
	3-0	2	0,45	ungefärbt/ undyed	PA 1156	24
	2-0	3	0,45	ungefärbt/ undyed	PA 1157	24
	4-0	1.5	0,70	ungefärbt/ undyed	PA 11517	24
	3-0	2	0,70	ungefärbt/ undyed	PA 11519	24
▼ DS 30 	0	3.5	0,70	violett/violet	PA 11411	24
	2-0	3	0,70	ungefärbt/ undyed	PA 1149	24
	3-0	2	0,70	ungefärbt/ undyed	PA 11421	24
	0	3.5	0,70	ungefärbt/ undyed	PA 1148	24
▼ DS 40 	2-0	3	0,70	ungefärbt/ undyed	PA 11420	24
1/2-Kreis gebogen außen schneidende = KS <i>1/2 curved, reverse cutting = KS</i>						
▼ KS 20 	3-0	2	0,45	violett/violet	PA 1181	24
Gerade schneidende = GS <i>straight, cutting = GS</i>						
▼ GS 12 	4-0	1.5	0,70	ungefärbt/ undyed	PA 1171	24
▼ GS 60 	2-0	3	0,70	ungefärbt/ undyed	PA 1172	24

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
1/2-Kreis schneidende Spatula, mit PREMIUM-Schliff = HSPM, HSPM F <i>1/4 circle, spatulated PREMIUM cut = HSPM, HSPM F</i>						
 2 x HSPM 5  PREMIUM-Schliff / PREMIUM cut	10-0	0.2	0,15	violett/violet	PA 1191	12
3/8-Kreis schneidende Spatula, mit PREMIUM-Schliff = DSPM <i>3/8 circle, spatulated, PREMIUM cut = DSPM</i>						
 DSPM 6  PREMIUM-Schliff / PREMIUM cut	9-0	0.3	0,30	violett/violet	PA 1203	12
	8-0	0.4	0,30	violett/violet	PA 1204	12
	7-0	0.5	0,45	violett/violet	PA 1205	24
 2 x DSPM 6  PREMIUM-Schliff / PREMIUM cut	9-0	0.3	0,30	violett/violet	PA 1206	12
	8-0	0.4	0,45	violett/violet	PA 1207	12
	7-0	0.5	0,45	violett/violet	PA 1208	24
 2 x DSPM 7  PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,45	violett/violet	PA 1201	24
 DSPM 9  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	violett/violet	PA 1202	24
1/4-Kreis schneidende Spatula, mit PREMIUM-Schliff = VSPM, VSPM F <i>1/4 circle, spatulated, PREMIUM cut = VSPM, VSPM F</i>						
 VSPM 6  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	violett/violet	PA 1214	24
 2 x VSPM 8  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	violett/violet	PA 12112	24
	6-0	0.7	0,45	violett/violet	PA 12113	24
 2 x VSPM 8 F  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	violett/violet	PA 1210	24
	5-0	1	0,45	violett/violet	PA 1211	24
	6-0	0.7	0,45	ungefärbt/ undyed	PA 1212	24
	5-0	1	0,45	ungefärbt/ undyed	PA 1213	24

PGA RESORBA®

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-needed

Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
USP	metric				
6-0	0.7	3 x 0,45	violett/violet	PS 133	24
5-0	1	3 x 0,45	violett/violet	PS 233	24
4-0	1.5	3 x 0,45	violett/violet	PS 333	24
3-0	2	3 x 0,45	violett/violet	PS 433	24
2-0	3	3 x 0,45	violett/violet	PS 533	24
0	3.5	3 x 0,45	violett/violet	PS 633	24
1	4	3 x 0,45	violett/violet	PS 733	24
2	5	3 x 0,45	violett/violet	PS 833	24
5	7	3 x 0,45	violett/violet	PS 1033	24
5-0	1	6 x 0,45	violett/violet	PS 265	24
4-0	1.5	6 x 0,45	violett/violet	PS 365	24
3-0	2	6 x 0,45	violett/violet	PS 465	24
2-0	3	6 x 0,45	violett/violet	PS 565	24
0	3.5	6 x 0,45	violett/violet	PS 665	24
1	4	6 x 0,45	violett/violet	PS 765	24
2	5	6 x 0,45	violett/violet	PS 865	24
4-0	1.5	12 x 0,45	violett/violet	PS 3125	24
3-0	2	12 x 0,45	violett/violet	PS 4125	24
2-0	3	12 x 0,45	violett/violet	PS 5125	24
0	3.5	12 x 0,45	violett/violet	PS 6125	24
1	4	12 x 0,45	violett/violet	PS 7125	24
2	5	12 x 0,45	violett/violet	PS 8125	24
5-0	1	2 x 0,70	violett/violet	PS 220	24
4-0	1.5	2 x 0,70	violett/violet	PS 320	24
3-0	2	2 x 0,70	violett/violet	PS 420	24
2-0	3	2 x 0,70	violett/violet	PS 520	24
0	3.5	2 x 0,70	violett/violet	PS 620	24
1	4	2 x 0,70	violett/violet	PS 720	24
2	5	2 x 0,70	violett/violet	PS 820	24
4-0	1.5	5 x 0,70	violett/violet	PS 350	24
3-0	2	5 x 0,70	violett/violet	PS 450	24
2-0	3	5 x 0,70	violett/violet	PS 550	24
0	3.5	5 x 0,70	violett/violet	PS 650	24
1	4	5 x 0,70	violett/violet	PS 750	24
2	5	5 x 0,70	violett/violet	PS 850	24
4	6	4 x 0,90	violett/violet	PS 949	24
5	7	4 x 0,90	violett/violet	PS 1049	24
4	6	3 x 1,00	violett/violet	PS 939	12
5	7	3 x 1,00	violett/violet	PS 1039	12
4-0	1.5	1 x 1,50	violett/violet	PS 315	24
3-0	2	1 x 1,50	violett/violet	PS 415	24
2-0	3	1 x 1,50	violett/violet	PS 515	24
0	3.5	1 x 1,50	violett/violet	PS 615	24
1	4	1 x 1,50	violett/violet	PS 715	24
2	5	1 x 1,50	violett/violet	PS 815	24
4	6	1 x 1,50	violett/violet	PS 915	24
5	7	1 x 1,50	violett/violet	PS 1015	24

PGA RESORBA®

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-needled

Stärke Nr. gauge no. USP metric		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
5-0	1	3 x 0,45	ungefärbt/ undyed	PU 233	24
4-0	1.5	3 x 0,45	ungefärbt/ undyed	PU 333	24
3-0	2	3 x 0,45	ungefärbt/ undyed	PU 433	24
2-0	3	3 x 0,45	ungefärbt/ undyed	PU 533	24
0	3.5	3 x 0,45	ungefärbt/ undyed	PU 633	24
1	4	3 x 0,45	ungefärbt/ undyed	PU 733	24
2	5	3 x 0,45	ungefärbt/ undyed	PU 833	24
4-0	1.5	6 x 0,45	ungefärbt/ undyed	PU 365	24
3-0	2	6 x 0,45	ungefärbt/ undyed	PU 465	24
2-0	3	6 x 0,45	ungefärbt/ undyed	PU 565	24
0	3.5	6 x 0,45	ungefärbt/ undyed	PU 665	24
1	4	6 x 0,45	ungefärbt/ undyed	PU 765	24
2	5	6 x 0,45	ungefärbt/ undyed	PU 865	24
3-0	2	12 x 0,45	ungefärbt/ undyed	PU 4125	24
2-0	3	12 x 0,45	ungefärbt/ undyed	PU 5125	24
0	3.5	12 x 0,45	ungefärbt/ undyed	PU 6125	24
1	4	12 x 0,45	ungefärbt/ undyed	PU 7125	24
2	5	12 x 0,45	ungefärbt/ undyed	PU 8125	24
2-0	3	2 x 0,70	ungefärbt/ undyed	PU 520	24
0	3.5	2 x 0,70	ungefärbt/ undyed	PU 620	24
2	5	2 x 0,70	ungefärbt/ undyed	PU 820	24
3-0	2	5 x 0,70	ungefärbt/ undyed	PU 450	24
2-0	3	5 x 0,70	ungefärbt/ undyed	PU 550	24
0	3.5	5 x 0,70	ungefärbt/ undyed	PU 650	24
1	4	5 x 0,70	ungefärbt/ undyed	PU 750	24
2	5	5 x 0,70	ungefärbt/ undyed	PU 850	24

PGA RESORBA®

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-needed

Stärke Nr. gauge no. USP metric		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
4-0	1.5	1 x 2,50	violett/violet	PL 100	12
3-0	2	1 x 2,50	violett/violet	PL 101	12
2-0	3	1 x 2,50	violett/violet	PL 102	12
0	3.5	1 x 2,50	violett/violet	PL 103	12
1	4	1 x 2,50	violett/violet	PL 104	12
2	5	1 x 2,50	violett/violet	PL 105	12
4	6	1 x 2,50	violett/violet	PL 106	12
5	7	1 x 2,50	violett/violet	PL 107	12

Sterile PGA-Schlingen, einzeln verpackt

Sterile PGA-Loops, single packed

Stärke Nr. gauge no. USP metric		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content
2-0	3	0,60	violett/violet	PSV 5	24
0	3.5	0,60	violett/violet	PSV 6	24
1	4	0,60	violett/violet	PSV 7	24
2	5	0,60	violett/violet	PSV 8	24