

RESOLON®



RESOLON® ist ein nicht-resorbierbares, monofiles Nahtmaterial aus Polyamid. Spezielle chemische Komponenten sorgen für geschmeidige Eigenschaften bei hoher Reißfestigkeit, was sich besonders bei der Hautnaht als vorteilhaft erwiesen hat.

- nicht-resorbierbares Material
- überdurchschnittlich geschmeidig und weich
- sehr gute Knüpfenschaften
- hervorragender Gewebedurchzug
- keine Kapillarität
- geringe Gewebereaktion



RESOLON®












RESOLON® is a non-absorbable monofilament suture based on polyamid. Special chemical components deliver smooth suture properties and high tensile strength supporting efficient skin closure.

- *non-absorbable suture*
- *softer and more supple than average*
- *perfect handling and knotting*
- *easy passage through tissue*
- *no capillarity*
- *non-irritating*

RESOLON®

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
1/2-Kreis schneidende Rundkörper = HRT <i>1/2 circle round bodied, trocar point = HRT</i>						
⊙ HRT 18 	5-0	1	0,45	blau/blue	88062	36
	4-0	1.5	0,45	blau/blue	88061	36
1/2-Kreis schneidende = HS, mit PREMIUM-Schliff = HSM <i>1/2 circle, reverse cutting = HS, PREMIUM cut = HSM</i>						
▼ HS 12 	5-0	1	0,45	blau/blue	88111	36
▽ HSM 17  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	blau/blue	88121	36
▼ HS 18 	5-0	1	0,45	blau/blue	88112	36
	4-0	1.5	0,45	blau/blue	88113	36
	4-0	1.5	0,75	blau/blue	88115	36
▼ HS 22 	3-0	2	0,45	blau/blue	88114	36
3/8-Kreis schneidende = DS, mit PREMIUM-Schliff = DSM, DSM F <i>3/8 circle, reverse cutting = DS, PREMIUM cut = DSM, DSM F</i>						
▽ DSM 7  PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,45	blau/blue	881543	36
▽ DSM 11  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	blau/blue	88151	36
	5-0	1	0,45	blau/blue	88152	36
▽ DSM 11 F  PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,45	blau/blue	881528	36
	6-0	0.7	0,45	blau/blue	881549	36
▼ DS 13 	6-0	0.7	0,45	blau/blue	881430	36
	5-0	1	0,45	blau/blue	881431	36
	4-0	1.5	0,45	blau/blue	881432	36

RESOLON®

Sterile Nadel-Faden-Kombination






Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 DSM 13  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	blau/blue	88153	36
	5-0	1	0,45	blau/blue	88154	36
	4-0	1.5	0,45	blau/blue	881523	36
	6-0	0.7	0,75	blau/blue	881553	36
	5-0	1	0,75	blau/blue	881554	36
 DS 16	6-0	0.7	0,45	blau/blue	881440	36
	5-0	1	0,45	blau/blue	88141	36
	4-0	1.5	0,45	blau/blue	88142	36
	3-0	2	0,45	blau/blue	881433	36
	5-0	1	0,75	blau/blue	881451	36
	4-0	1.5	0,75	blau/blue	881452	36
	6-0	0.7	0,45	blau/blue	881529	36
 DSM 16 PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	blau/blue	88155	36
	4-0	1.5	0,45	blau/blue	88156	36
	3-0	2	0,45	blau/blue	881548	36
	5-0	1	0,75	blau/blue	881530	36
	4-0	1.5	0,75	blau/blue	881551	36
	3-0	2	0,75	blau/blue	881552	36
	5-0	1	0,45	blau/blue	881550	36
 DSM 16 F PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	blau/blue	881550	36
 DS 18	5-0	1	0,45	blau/blue	88143	36
	4-0	1.5	0,45	blau/blue	88144	36
	3-0	2	0,45	blau/blue	88146	36
	2-0	3	0,45	blau/blue	881454	36
	5-0	1	0,75	blau/blue	881450	36
	4-0	1.5	0,75	blau/blue	88145	36
	3-0	2	0,75	blau/blue	88147	36
	2-0	3	0,75	blau/blue	88148	36
	4-0	1.5	4 x 0,45	blau/blue	M 881421 Multipack	24
	4-0	1.5	0,45	blau/blue	881441	36
 DS 18 F	4-0	1.5	0,45	blau/blue	881441	36
	4-0	1.5	0,75	blau/blue	881442	36

RESOLON®

Sterile Nadel-Faden-Kombination






Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 DSM 18 PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	blau/blue	88157	36
	4-0	1.5	0,45	blau/blue	88158	36
	3-0	2	0,45	blau/blue	881510	36
	2-0	3	0,45	blau/blue	881531	36
	5-0	1	0,75	blau/blue	881555	36
	4-0	1.5	0,75	blau/blue	88159	36
	3-0	2	0,75	blau/blue	881511	36
	2-0	3	0,75	blau/blue	881512	36
	 DSM 18 F PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	blau/blue	881541
4-0		1.5	0,75	blau/blue	881542	36
4-0		1.5	4 x 0,45	blau/blue	M 881548	24
					Multipack	
 DS 21	4-0	1.5	0,45	blau/blue	88149	36
	3-0	2	0,45	blau/blue	881411	36
	4-0	1.5	0,75	blau/blue	881410	36
	3-0	2	0,75	blau/blue	881412	36
	2-0	3	0,75	blau/blue	881413	36
 DSM 21 PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	blau/blue	881513	36
	3-0	2	0,45	blau/blue	881515	36
	4-0	1.5	0,75	blau/blue	881514	36
	3-0	2	0,75	blau/blue	881516	36
	2-0	3	0,75	blau/blue	881517	36
 DS 24	5-0	1	0,45	blau/blue	881458	36
	4-0	1.5	0,45	blau/blue	881424	36
	3-0	2	0,45	blau/blue	881414	36
	2-0	3	0,45	blau/blue	881425	36
	4-0	1.5	0,75	blau/blue	881457	36
	3-0	2	0,75	blau/blue	881415	36
	2-0	3	0,75	blau/blue	881416	36

RESOLON®

Sterile Nadel-Faden-Kombination






Needle - suture - combination, sterile

nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 DSM 24 PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	blau/blue	881544	36
	3-0	2	0,45	blau/blue	881518	36
	2-0	3	0,45	blau/blue	881524	36
	4-0	1.5	0,75	blau/blue	881545	36
	3-0	2	0,75	blau/blue	881519	36
	2-0	3	0,75	blau/blue	881520	36
 DS 27	3-0	2	0,45	blau/blue	881453	36
	2-0	3	0,45	blau/blue	881423	36
	0	3.5	0,45	blau/blue	881456	36
	3-0	2	0,75	blau/blue	881417	36
	2-0	3	0,75	blau/blue	881418	36
	0	3.5	0,75	blau/blue	881556	36
 DSM 27 PREMIUM-Schliff / PREMIUM cut	2-0	3	0,45	blau/blue	881525	36
	2-0	3	0,75	blau/blue	881526	36
 DS 30	4-0	1.5	0,45	blau/blue	881434	36
	3-0	2	0,45	blau/blue	881426	36
	2-0	3	0,45	blau/blue	881428	36
	0	3.5	0,45	blau/blue	881429	36
	4-0	1.5	0,75	blau/blue	881435	36
	3-0	2	0,75	blau/blue	881427	36
	2-0	3	0,75	blau/blue	881419	36
	0	3.5	0,75	blau/blue	881420	36
 DSM 30 PREMIUM-Schliff / PREMIUM cut	3-0	2	0,45	blau/blue	881532	36
	2-0	3	0,45	blau/blue	881534	36
	0	3.5	0,45	blau/blue	881536	36
	3-0	2	0,75	blau/blue	881533	36
	2-0	3	0,75	blau/blue	881535	36
	0	3.5	0,75	blau/blue	881537	36

RESOLON®

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ DS 40 	2-0	3	0,45	blau/blue	881436	36
	0	3.5	0,45	blau/blue	881438	36
	2-0	3	0,75	blau/blue	881437	36
	0	3.5	0,75	blau/blue	881439	36
	2-0	3	0,90	blau/blue	881455	36
▽ DSM 40  PREMIUM-Schliff / PREMIUM cut	2-0	3	0,75	blau/blue	881538	36
	0	3.5	0,75	blau/blue	881540	36
Gerade Rundkörper = GRW <i>Straight round bodied = GRW</i>						
● GR 20 	4-0	1.5	0,45	blau/blue	88093	36
Gerade Rundkörper, biegsame Nadel = GRW <i>Straight round bodied, flexible needle = GRW</i>						
● 2 x GRW 65  Biegsame Nadel / Flexible needle	2-0	3	0,75	blau/blue	88091	36
	0	3.5	0,75	blau/blue	88092	36
Gerade schneidende = GS <i>straight, cutting = GS</i>						
▼ GS 60 	3-0	2	0,75	blau/blue	88171	36
	2-0	3	0,75	blau/blue	88172	36

RESOLON®

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-needed

Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
USP	metric				
4-0	1.5	3 x 0,50	blau/blue	X 335 N	24
3-0	2	3 x 0,50	blau/blue	X 435 N	24
2-0	3	3 x 0,50	blau/blue	X 535 N	24
0	3.5	3 x 0,50	blau/blue	X 635 N	24
3-0	2	10 x 0,50	blau/blue	X 4105 N	24