

SEIDE



SEIDE ist ein natürliches, nicht-resorbierbares Nahtmaterial, welches aus dem Kokon der Seidenspinnerraupe gewonnen wird. Die nachfolgende Veredlung sorgt für gute Handlungseigenschaften und bietet vielfältige Einsatzmöglichkeiten. SEIDE wird in schwarz und weiß angeboten sowie als gezwirnte Variante in blau.

- nicht-resorbierbarer Seidenfaden
- hohe Geschmeidigkeit
- hervorragendes Knüpfverhalten
- maximale Zugfestigkeit



SILK











SEIDE is a natural non-absorbable suture that is made from the continuous fibres from silkworm cocoons. After purifying the raw material the product offers good handling properties and multiple indications to use. SEIDE is available in black or white and also blue in a twisted version.

- *non-absorbable suture*
- *permanent pliability*
- *perfect knot security*
- *maximum tensile strength*

SEIDE

Sterile Nadel-Faden-Kombination







Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
1/2-Kreis Rundkörper = HR, HRS <i>1/2 circle, round bodied = HR, HRS</i>						
● HR 10 	5-0	1	0,45	schwarz, geflochten <i>black, braided</i>	4021	36
● HR 12 	5-0	1	0,45	schwarz, geflochten <i>black, braided</i>	4023	36
	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	4024	36
● HR 12 A  Abziehnadel/ <i>control release needle</i>	4-0	1.5	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40248 Multipack	12
	3-0	2	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40264 Multipack	12
● HR 17 	5-0	1	0,75	schwarz, geflochten <i>black, braided</i>	4026	36
	4-0	1.5	0,75	schwarz, geflochten <i>black, braided</i>	4027	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	40267	36
	3-0	2	0,75	schwarz, geflochten <i>black, braided</i>	4028	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	4029	36
● HR 17 A  Abziehnadel/ <i>control release needle</i>	4-0	1.5	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40249 Multipack	12
	3-0	2	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40250 Multipack	12
	2-0	3	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40265 Multipack	12
● 2 x HR 17 	4-0	1.5	0,75	schwarz, geflochten <i>black, braided</i>	40215	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	40216	36
● HR 22 	4-0	1.5	0,75	schwarz, geflochten <i>black, braided</i>	40219	36
	3-0	2	0,75	schwarz, geflochten <i>black, braided</i>	40220	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	40221	36
● HR 22 A  Abziehnadel/ <i>control release needle</i>	3-0	2	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40252 Multipack	12
	2-0	3	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A 40253 Multipack	12

SEIDE

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● HR 27 	4-0	1.5	0,75	schwarz, geflochten <i>black, braided</i>	40225	36
	3-0	2	0,75	schwarz, geflochten <i>black, braided</i>	40226	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	40227	36
	0	3.5	0,75	schwarz, geflochten <i>black, braided</i>	40228	36
	1	4	0,75	schwarz, geflochten <i>black, braided</i>	40229	36
● HR 27 A  Abziehnadel/ <i>control release needle</i>	2-0	3	8 x 0,45	schwarz, geflochten <i>black, braided</i>	A40259 Multipack	12
● HR 30 	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	40238	36
	0	3.5	0,75	schwarz, geflochten <i>black, braided</i>	40239	36
● HR 35 	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	40242	36
	0	3.5	0,75	schwarz, geflochten <i>black, braided</i>	40243	36
	1	4	0,75	schwarz, geflochten <i>black, braided</i>	40244	36
	2	5	0,75	schwarz, geflochten <i>black, braided</i>	40245	36
● HRS 35 	0	3.5	0,75	schwarz, geflochten <i>black, braided</i>	4031	36
● HR 40 	1	4	0,75	schwarz, geflochten <i>black, braided</i>	40247	36

SEIDE

Sterile Nadel-Faden-Kombination










Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
1/2-Kreis schneidende Rundkörper = HRT <i>1/2 circle, round bodied, trocar point = HRT</i>						
▼ HRT 18 	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	40616	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	40615	36
	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	40617	36
	3-0	2	6 x 0,45	schwarz, geflochten <i>black, braided</i>	M 40618 Multipack	12
	4-0	1.5	0,75	schwarz, geflochten <i>black, braided</i>	4062	36
	3-0	2	0,75	schwarz, geflochten <i>black, braided</i>	4063	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	4064	36
	▼ HRT 25 	3-0	2	0,75	schwarz, geflochten <i>black, braided</i>	40610
2-0		3	0,75	schwarz, geflochten <i>black, braided</i>	40611	36
Asymptotische schneidende Rundkörper = ART <i>Asymptotic, round bodied, trocar point = ART</i>						
▼ ART 13 	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	4281	12
	▼ ART 25 	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	4282
3-0		2	0,45	schwarz, geflochten <i>black, braided</i>	4283	12
3/8-Kreis Rundkörper = DR, DR F <i>3/8 circle, round bodied = DR, DR F</i>						
● DR 10 F 	7-0	0.5	0,30	schwarz, geflochten <i>black, braided</i>	4071	24
	6-0	0.7	0,45	schwarz, geflochten <i>black, braided</i>	4072	24
● DR 12 	5-0	1	0,75	schwarz, geflochten <i>black, braided</i>	4078	36
	4-0	1.5	0,75	schwarz, geflochten <i>black, braided</i>	4079	36
● 2 x DR 12 	5-0	1	0,45	schwarz, geflochten <i>black, braided</i>	40712	36

SEIDE

Sterile Nadel-Faden-Kombination










Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
● DR 12 F 	6-0	0.7	0,45	schwarz, geflochten black, braided	4077	36
● DR 18 	3-0	2	0,30	schwarz, geflochten black, braided	40725	36
3/8-Kreis Rundkörper stumpf 3/8 circle, round bodied, blunt point						
○ DRN 17 	3-0	2	0,30	schwarz, geflochten black, braided	4273	36
	3-0	2	2 x 0,75	schwarz, geflochten black, braided	M 4272 Multipack	24
3/8-Kreis schneidende Rundkörper = DRT 3/8 circle, round bodied, trocar point = DRT						
⊖ DRT 18 	3-0	2	0,75	schwarz, geflochten black, braided	4083	36
Gerade Rundkörper = GR Straight round bodied = GR						
● GR 20 	4-0	1.5	0,45	schwarz, geflochten black, braided	40918	36
	4-0	1.5	0,75	schwarz, geflochten black, braided	4091	36
● GR 45 	3-0	2	0,75	schwarz, geflochten black, braided	4095	36
1/2-Kreis schneidende = HS mit PREMIUM-Schliff = HSM, HSM S 1/2 circle, reverse cutting = HS PREMIUM cut = HSM, HSM S						
▽ 2 x HSM 9 S  PREMIUM-Schliff PREMIUM cut	4-0	1.5	0,45	schwarz, geflochten black, braided	4125	24
▽ 2 x HSM 10  PREMIUM-Schliff PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten black, braided	4130	24
▼ HS 12 	5-0	1	0,45	schwarz, geflochten black, braided	4111	36

SEIDE

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▽ HSM 13  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	schwarz, geflochten black, braided	41210	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	41211	36
▼ HS 18 	5-0	1	0,45	schwarz, geflochten black, braided	4118	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	4112	36
	3-0	2	0,45	schwarz, geflochten black, braided	4113	36
	4-0	1.5	0,75	schwarz, geflochten black, braided	4119	36
	3-0	2	0,75	schwarz, geflochten black, braided	41111	36
▼ HS 22 	3-0	2	0,45	schwarz, geflochten black, braided	4116	36
	2-0	3	0,45	schwarz, geflochten black, braided	4117	36
▼ HS 35 	2-0	3	0,75	schwarz, geflochten black, braided	41110	36
3/8-Kreis schneidende = DS mit PREMIUM-Schliff = DSM, DSM F <i>3/8 circle, reverse cutting = DS PREMIUM cut = DSM, DSM F</i>						
▽ DSM 6  PREMIUM-Schliff / PREMIUM cut	8-0	0.4	0,30	blau, gedreht blue, twisted	4154	24
▽ DSM 7 F  PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,30	schwarz, geflochten black, braided	4157	24
	6-0	0.7	0,30	schwarz, geflochten black, braided	4158	24
▽ 2 x DSM 7 F  PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,45	schwarz, geflochten black, braided	4159	24
	6-0	0.7	0,45	schwarz, geflochten black, braided	41510	24
▽ DSM 11  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten black, braided	41517	36
▽ DSM 11 F  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten black, braided	41538	36

SEIDE

Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 2 x DSM 11 F PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten black, braided	41519	36
 DS 13	6-0	0.7	0,35	schwarz, geflochten black, braided	4141	36
	5-0	1	0,45	schwarz, geflochten black, braided	4142	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	4143	36
 2 x DS 13	5-0	1	0,35	schwarz, geflochten black, braided	4145	36
	4-0	1.5	0,35	schwarz, geflochten black, braided	4146	36
	6-0	0.7	0,45	schwarz, geflochten black, braided	4144	36
 DSM 13 PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten black, braided	41542	36
	5-0	1	0,45	schwarz, geflochten black, braided	41520	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	41521	36
 DS 16	6-0	0.7	0,45	schwarz, geflochten black, braided	4147	36
	5-0	1	0,45	schwarz, geflochten black, braided	4148	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	4149	36
 DSM 16 PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	schwarz, geflochten black, braided	41524	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	41525	36
 DS 18	5-0	1	0,45	schwarz, geflochten black, braided	41411	36
	4-0	1.5	0,45	schwarz, geflochten black, braided	41412	36
	3-0	2	0,45	schwarz, geflochten black, braided	41413	36
	4-0	1.5	0,75	schwarz, geflochten black, braided	41414	36
	3-0	2	0,75	schwarz, geflochten black, braided	41438	36

SEIDE

Sterile Nadel-Faden-Kombination














Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
 DSM 18 PREMIUM-Schliff / <i>PREMIUM cut</i>	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	41527	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	41528	36
	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41529	36
 DSM 18 F PREMIUM-Schliff / <i>PREMIUM cut</i>	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	41540	36
 DS 21	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	41431	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	41432	36
	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41433	36
 DSM 21 PREMIUM-Schliff / <i>PREMIUM cut</i>	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	41534	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	41535	36
	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41536	36
 DS 24	4-0	1.5	0,45	schwarz, geflochten <i>black, braided</i>	41416	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	41417	36
	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41418	36
	3-0	2	0,75	schwarz, geflochten <i>black, braided</i>	41419	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	41436	36
 DSM 24 PREMIUM-Schliff / <i>PREMIUM cut</i>	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41533	36
	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	41532	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	41531	36
 DS 27	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41422	36
 DS 30	3-0	2	0,45	schwarz, geflochten <i>black, braided</i>	41423	36
	2-0	3	0,45	schwarz, geflochten <i>black, braided</i>	41424	36
	0	3.5	0,45	schwarz, geflochten <i>black, braided</i>	41425	36
	2-0	3	0,75	schwarz, geflochten <i>black, braided</i>	41426	36
	0	3.5	0,75	schwarz, geflochten <i>black, braided</i>	41428	36
	1	4	0,75	schwarz, geflochten <i>black, braided</i>	41439	36

SEIDE

Sterile Nadel-Faden-Kombination







Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ DS 40 	0	3.5	0,75	schwarz, geflochten black, braided	41430	36
	1	4	0,75	schwarz, geflochten black, braided	41437	36
Gerade schneidende = GS <i>Straight, cutting = GS</i>						
▼ GS 60 	2-0	3	0,75	schwarz, geflochten black, braided	4170	36
1/2-Kreis gebogene schneidende = KS <i>1/2 circle, curved reverse cutting = KS</i>						
▼ KS 20 	3-0	2	0,45	schwarz, geflochten black, braided	4181	36
1/2-Kreis schneidende Spatula, PREMIUM-Schliff = HSPM, HSPM F <i>1/2 circle, spatulated, PREMIUM cut = HSPM, HSPM F</i>						
 2 x HSPM 6 F  PREMIUM-Schliff PREMIUM cut	8-0	0.4	0,30	blau, gedreht blue, twisted	4195	24
 2 x HSPM 8  PREMIUM-Schliff PREMIUM cut	6-0	0.7	0,45	weiß, geflochten white, braided	4191	24
3/8-Kreis schneidende Spatula, PREMIUM-Schliff = DSPM <i>3/8 circle, spatulated, PREMIUM cut = DSPM</i>						
 DSPM 5  PREMIUM-Schliff PREMIUM cut	8-0	0.4	0,30	blau, gedreht blue, twisted	4202	24
 2 x DSPM 5  PREMIUM-Schliff PREMIUM cut	8-0	0.4	0,30	blau, gedreht blue, twisted	4206	24
 2 x DSPM 6  PREMIUM-Schliff PREMIUM cut	8-0	0.4	0,30	blau, gedreht blue, twisted	42011	24
	9-0	0.3	0,45	blau, gedreht blue, twisted	42017	24

SEIDE

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
  PREMIUM-Schliff / PREMIUM cut	8-0	0.4	0,45	schwarz, geflochten <i>black, braided</i>	42015	24
	8-0	0.4	0,45	blau, gedreht <i>black, braided</i>	42024	24
	7-0	0.5	0,45	schwarz, geflochten <i>black, braided</i>	42021	24
  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten <i>black, braided</i>	4208	24
1/4-Kreis schneidende Spatula, PREMIUM-Schliff = VSPM F <i>1/4 circle, spatulated, PREMIUM cut = VSPM F</i>						
  PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz, geflochten <i>black, braided</i>	4217	24
	5-0	1	0,45	schwarz, geflochten <i>black, braided</i>	4218	24

SEIDE

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-neededled

Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
USP	metric				
4-0	1.5	5 x 0,35	schwarz/black	G 353	24
3-0	2	5 x 0,35	schwarz/black	G 453	24
2-0	3	5 x 0,35	schwarz/black	G 553	24
4-0	1.5	3 x 0,50	schwarz/black	G 335	24
3-0	2	3 x 0,50	schwarz/black	G 435	24
2-0	3	3 x 0,50	schwarz/black	G 535	24
0	3.5	3 x 0,50	schwarz/black	G 635	24
1	4	3 x 0,50	schwarz/black	G 735	24
2	5	3 x 0,50	schwarz/black	G 835	24
0	3.5	3 x 0,50	weiß/white	WG 635	24
1	4	3 x 0,50	weiß/white	WG 735	24
4-0	1.5	5 x 0,50	schwarz/black	G 355	24
3-0	2	5 x 0,50	schwarz/black	G 455	24
2-0	3	5 x 0,50	schwarz/black	G 555	24
0	3.5	5 x 0,50	schwarz/black	G 655	24
1	4	5 x 0,50	schwarz/black	G 755	24
2	5	5 x 0,50	schwarz/black	G 855	24
3	6	5 x 0,50	schwarz/black	G 955	24
2-0	3	5 x 0,50	weiß/white	WG 555	24
0	3.5	5 x 0,50	weiß/white	WG 655	24
3-0	2	6 x 0,50	schwarz/black	SG 455	24
2-0	3	6 x 0,50	schwarz/black	SG 555	24
0	3.5	6 x 0,50	schwarz/black	SG 655	24
1	4	6 x 0,50	schwarz/black	SG 755	24
2	5	6 x 0,50	schwarz/black	SG 855	24
2	5	6 x 0,50	weiß/white	WSG 855	24
5-0	1	10 x 0,50	schwarz/black	G 2105	24
4-0	1.5	10 x 0,50	schwarz/black	G 3105	24
3-0	2	10 x 0,50	schwarz/black	G 4105	24
2-0	3	10 x 0,50	schwarz/black	G 5105	24
0	3.5	10 x 0,50	schwarz/black	G 6105	24
1	4	10 x 0,50	schwarz/black	G 7105	24
2-0	3	2 x 0,75	schwarz/black	G 520	24
0	3.5	2 x 0,75	schwarz/black	G 620	24
1	4	2 x 0,75	schwarz/black	G 720	24
2	5	2 x 0,75	schwarz/black	G 820	24
3	6	2 x 0,75	schwarz/black	G 920	24
5	7	2 x 0,75	schwarz/black	G 1020	24
1	4	2 x 0,75	weiß/white	WG 720	24
2	5	2 x 0,75	weiß/white	WG 820	24
3	6	2 x 0,75	weiß/white	WG 920	24
3-0	2	5 x 0,75	schwarz/black	G 450	24
2-0	3	5 x 0,75	schwarz/black	G 550	24
0	3.5	5 x 0,75	schwarz/black	G 650	24
1	4	5 x 0,75	schwarz/black	G 750	24
2	5	5 x 0,75	schwarz/black	G 850	24
3-0	2	1 x 1,40	schwarz/black	SG 411	24
2-0	3	1 x 1,40	schwarz/black	SG 511	24
5	7	1 x 1,40	schwarz/black	SG 1011	24